# Tcf4 (TCF7L2) Antibody Staining Embryonic muscle cryosections Kardon Lab April 2012 (Allyson Merrell)

### **Protocol:**

Embed embryos in OCT and cryosection (8-12 um thick)

#### **Day 1:**

5 min 4% PFA/PBS 3x3 min PBS 1 hour block in 5% goat or donkey serum/PBS at room temperature Overnight at 4°C in 1° antibody

## **Day 2:**

5x3 min PBS
2 hours 2<sup>0</sup> antibody at room temperature
3x3 min PBS
1 hour in Strep at room temperature
3x3 min PBS
5 min 4% PFA/PBS
3x3 min PBS
5 min Hoechst
2x3 min PBS
Rinse H<sub>2</sub>O and coverslip with fluoromount

### Reagents:

Tcf-4 (TCF7L2) Antibodies:

Rabbit mAb (1:100=0.7ug/ml final concentration): Clone C48H11 Cell Signaling Technology 2569

Mouse IgG2a (1:100=10ug/ml final concentration): Clone 6H5-3 Millipore 05-511 Rabbit antibody is preferred

Biotin conjugated goat α rabbit

Jackson ImmunoResearch 111-065-144, stock diluted 1:1 in glycerol 1:500 (1.2  $\mu$ g/ml) in 5% goat serum in PBS

Strep: Streptavidin conjugated Dylight 594 or 488 (Jackson ImmunoResearch Dylight 594: 016-510-084, Dylight 488: 016-480-084, stock diluted 1:1 in glycerol)

1:200 (4.5 µg/ml) in PBS

Hoechst (Invitrogen, H1399) 10ul in 250 ml PBS

Fluoromount-G (Southern Biotech #0100-01)

All PBS washes and Hoechst staining in 250 ml buckets

 $1^{\circ}$  and  $2^{\circ}$  antibodies and block in small volume (200-500ul) on slide, laying flat in humidified darkened plastic tray